



BacGO

P029
1 μmol

- **Known Property** gram positive bacteria probe
- **Application** Immunofluorescence
- **Cell selectivity mechanism:** COLD (peptidoglycan of bacteria)
- **Storage**
 - ① Delivery: Room Temperature
 - ② Dried compound: 4 °C or -20 °C
 - ③ Compound solution: 4 °C or -20 °C

ORDER

- SenPro
- order@senprobe.com
- www.senprobe.com


General Use Guide

More than 1/100 dilution of 10mM of DMSO stock solution is essential

For biomedical use to avoid DMSO concentration higher than 1%.

Working concentrations for specific applications should be determined by the investigator.

It is recommended to use up the buffer diluted solution within one day. The compound may be decomposed or precipitated out from buffer solution.




Molecular Weight

485.1 (C₂₆H₂₃B₂F₂N₃O₃)

$\lambda_{\text{ex}} / \lambda_{\text{em}}$

563 / 578 nm

BacGO (Bacteria Gram-positive Orange) is selective fluorescent probe for gram positive bacteria over gram negative bacteria. Boronic acid is a popular binding motif for vicinal diol and has been used for carbohydrate recognition. Gram positive bacteria contains higher contents of peptidoglycan than Gram negative bacteria and boronic acid is expected to more selectively bind to Gram positive bacteria. A series of boronic acid containing fluorescent compounds were tested for Gram positive bacteria selective staining, and **BacGO** was discovered as a universal Gram positive bacteria probe.



Universal selectivity of **BacGO** towards Gram-positive bacterial strains. (a) **BacGO** was incubated with 16 bacterial strains for 5 or 60 min. Flow-cytometry histogram images of **BacGO** 1 μm staining. (b) The structure of **BacGO**. (c) Images of Gram staining and fluorescent staining using 16 bacterial strains. The fluorescent image of **BacGO** (Texas-red filter) and Gram-staining images were taken using a Zeiss Microscope equipped with a 100 x objective lens. Scalebar = 10 μm . (BC, *Bacillus cereus*; BM, *Bacillus megaterium*; BS, *Bacillus subtilis*; CS, *Cronobacter sakazakii*; EC, *Escherichia coli*; ED, *Escherichia coli* DH5 α ; UTI, *Escherichia coli* UTI89; EF, *Enterococcus faecalis*; KP, *Klebsiella pneumoniae*; LF, *Lactobacillus fermentum*; LP, *Lactobacillus plantarum*; LS, *Lactobacillus sakei*; PA, *Pseudomonas aeruginosa*; PP, *Pseudomonas putida*; SA, *Staphylococcus aureus*; ST, *Streptococcus thermophilus*).

- Related probes: CDy11, CDy14, CDr15

Reference

1. **Development of a Universal Fluorescent Probe for Gram-Positive Bacteria**, Kwon, H. Y.; Liu, X.; Choi, E. G.; Lee, J. Y.; Choi, S. Y.; Kim, J. Y.; Wang, L.; Park, S. J.; Kim, B.; Lee, Y. A.; Kim, J. J.; Kang, N. Y.*; Chang, Y. T.* *Angew. Chem. Int. Ed. Engl.* 2019, 58, 8426-8431.